

SLC2A13

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:4000

Calculated MW: 70kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human SLC2A13

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #: A9993

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 114134

Isotype: IgG

Swiss Prot: Q96QE2

Purity: Affinity purification

For research use only.