

PELI3

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 51kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human PELI3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A9999

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 246330

Isotype: IgG

Swiss Prot: Q8N2H9

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a scaffold protein and an intermediate signaling protein in the innate immune response pathway. The encoded protein helps transmit the immune response signal from Toll-like receptors to IRAK1/TRAF6 complexes. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2011]

To place an order, please [Click HERE](#).