

ACTB

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 43kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ACTB

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

h

Synonym:

BRWS1; PS1TP5BP1;

Catalog #: AC006

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 60

Isotype: IgG

Swiss Prot: P60709

Purity: Affinity purification

For research use only.

Background:

This gene encodes one of six different actin proteins. Actins are highly conserved proteins that are involved in cell motility, structure, and integrity. This actin is a major constituent of the contractile apparatus and one of the two nonmuscle cytoskeletal actins.

To place an order, please [Click HERE](#).