

-Tubulin

Reactivity: Human Mouse Rat

Tested applications: WB IF

Observed MW: Refer to figures

Immunogen:

Recombinant Protein

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #: AC021

Antibody Type:

Monoclonal Antibody

Species: Mouse

Gene ID: 203068

Swiss Prot: P07437

For research use only.