

MYC-Tag

Tested applications:WB ICC IP FC

Recommended Dilution:WB 1:500 - 1:1000 ICC 1:50 - 1:200 IP 1:50 - 1:100 FC 1:50 - 1:100

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of EQKLISEEDL (MYC-tag)-KLH

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #:AE038

Antibody Type:

Monoclonal Antibody

Species:Mouse

Isotype:IgG

Purity:Affinity purification

For research use only.