

Phospho-HSP90AB1-S254

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 92kDa

Observed MW: Refer to Figures

Immunogen:

A phospho specific peptide corresponding to residues surrounding S254 of human HSP90AB1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

hj

Synonym:

HSP84; HSPC2; HSPCB; D6S182; HSP90B;

Catalog #: AP0203

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3326

Isotype: IgG

Swiss Prot: P08238

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the heat shock protein 90 family; these proteins are involved in signal transduction, protein folding and degradation and morphological evolution. This gene encodes the constitutive form of the cytosolic 90 kDa heat-shock protein and is thought to play a role in gastric apoptosis and inflammation. Alternative splicing results in multiple transcript variants. Pseudogenes have been identified on multiple chromosomes.

To place an order, please [Click HERE](#).