

## Phospho-MET-Y1349

**Reactivity:** Human Mouse Rat

**Tested applications:** WB

**Recommended Dilution:** WB 1:500 - 1:2000

**Calculated MW:** 156kDa

**Observed MW:** Refer to Figures

**Immunogen:**

A phospho specific peptide corresponding to residues surrounding Y1349 of human MET

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

HGFR; AUTS9; RCCP2; c-Met;

**Background:**

The proto-oncogene MET product is the hepatocyte growth factor receptor and encodes tyrosine-kinase activity. The primary single chain precursor protein is post-translationally cleaved to produce the alpha and beta subunits, which are disulfide linked to form the mature receptor. Various mutations in the MET gene are associated with papillary renal carcinoma. Two transcript variants encoding different isoforms have been found for this gene.

*To place an order, please [Click HERE](#).*

**Catalog #:** AP0230

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 4233

**Isotype:** IgG

**Swiss Prot:** P08581

**Purity:** Affinity purification

For research use only.