

Phospho-DOK1-Y362

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 62kDa

Observed MW: Refer to Figures

Immunogen:

A phospho specific peptide corresponding to residues surrounding Y362 of human DOK1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

P62DOK;

Catalog #: AP0248

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1796

Isotype: IgG

Swiss Prot: Q99704

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is part of a signal transduction pathway downstream of receptor tyrosine kinases. The encoded protein is a scaffold protein that helps form a platform for the assembly of multiprotein signaling complexes. Two transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).