

Phospho-ITGB3-Y785

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:1000

Calculated MW: 110kDa

Observed MW: Refer to Figures

Immunogen:

A phospho specific peptide corresponding to residues surrounding Y785 of human ITGB3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

q

Synonym:

GT; CD61; GP3A; BDPLT2; GPIIIa; BDPLT16;

Catalog #: AP0252

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3690

Isotype: IgG

Swiss Prot: P05106

Purity: Affinity purification

For research use only.

Background:

The ITGB3 protein product is the integrin beta chain beta 3. Integrins are integral cell-surface proteins composed of an alpha chain and a beta chain. A given chain may combine with multiple partners resulting in different integrins. Integrin beta 3 is found along with the alpha IIb chain in platelets. Integrins are known to participate in cell adhesion as well as cell-surface mediated signalling.

To place an order, please [Click HERE](#).