

Phospho-CTNNA1-S641

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 100kDa

Observed MW: Refer to Figures

Immunogen:

A phospho specific peptide corresponding to residues surrounding S641 of human CTNNA1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

jt

Synonym:

CAP102;

To place an order, please [Click HERE](#).

Catalog #: AP0273

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1495

Isotype: IgG

Swiss Prot: P35221

Purity: Affinity purification

For research use only.