

Phospho-EIF4B-S422

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 80kDa

Observed MW: Refer to Figures

Immunogen:

A phospho specific peptide corresponding to residues surrounding S422 of human EIF4B

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

jq

Synonym:

EIF-4B; PRO1843;

To place an order, please [Click HERE](#).

Catalog #: AP0278

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1975

Isotype: IgG

Swiss Prot: P23588

Purity: Affinity purification

For research use only.