

Phospho-AKT2-S474

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 60kDa

Observed MW: Refer to Figures

Immunogen:

A phospho specific peptide corresponding to residues surrounding S474 of human AKT2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

opr

Synonym:

PKBB; PRKBB; HIHGGH; PKBBETA; RAC-BETA;

Catalog #: AP0305

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 208

Isotype: IgG

Swiss Prot: P31751

Purity: Affinity purification

For research use only.

Background:

This gene is a putative oncogene encoding a protein belonging to a subfamily of serine/threonine kinases containing SH2-like (Src homology 2-like) domains. The gene was shown to be amplified and overexpressed in 2 of 8 ovarian carcinoma cell lines and 2 of 15 primary ovarian tumors.

Overexpression contributes to the malignant phenotype of a subset of human ductal pancreatic cancers. The encoded protein is a general protein kinase capable of phosphorylating several known proteins.

To place an order, please [Click HERE](#).