

Phospho-BCL2-T56

Reactivity: Human

Tested applications: IHC

Recommended Dilution: IHC 1:50 - 1:100

Calculated MW: 26kDa

Observed MW: Refer to Figures

Immunogen:

A phospho specific peptide corresponding to residues surrounding T56 of human BCL2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

bi

Synonym:

Bcl-2; PPP1R50;

Catalog #: AP0316

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 596

Isotype: IgG

Swiss Prot: P10415

Purity: Affinity purification

For research use only.

Background:

This gene encodes an integral outer mitochondrial membrane protein that blocks the apoptotic death of some cells such as lymphocytes. Constitutive expression of BCL2, such as in the case of translocation of BCL2 to Ig heavy chain locus, is thought to be the cause of follicular lymphoma.

Two transcript variants, produced by alternate splicing, differ in their C-terminal ends.

To place an order, please [Click HERE](#).