

Phospho-CFL1/CFL2-Y88

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:1000 IHC 1:50 - 1:100

Calculated MW: 19kDa

Observed MW: Refer to Figures

Immunogen:

A phospho specific peptide corresponding to residues surrounding Y88 of human CFL1/CFL2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CFL1/CFL2;

Catalog #: AP0329

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1072/ 1073

Isotype: IgG

Swiss Prot: P23528/Q9Y281

Purity: Affinity purification

For research use only.

Background:

Controls reversibly actin polymerization and depolymerization in a pH-sensitive manner. It has the ability to bind G- and F-actin in a 1:1 ratio of cofilin to actin. It is the major component of intranuclear and cytoplasmic actin rods.

To place an order, please [Click HERE](#).