

Phospho-FAM129A-S602

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:1000

Calculated MW: 150kDa

Observed MW: Refer to Figures

Immunogen:

A phospho specific peptide corresponding to residues surrounding S602 of human FAM129A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

FAM129A;

Catalog #: AP0350

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 116496

Isotype: IgG

Swiss Prot: Q9BZQ8

Purity: Affinity purification

For research use only.

Background:

Regulates phosphorylation of a number of proteins involved in translation regulation including EIF2A, EIF4EBP1 and RPS6KB1. May be involved in the endoplasmic reticulum stress response

To place an order, please [Click HERE](#).