

Phospho-HSF1-S303

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 82kDa

Observed MW: Refer to Figures

Immunogen:

A phospho specific peptide corresponding to residues surrounding S303 of human HSF1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

k

Synonym:

HSTF1;

Catalog #: AP0363

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3297

Isotype: IgG

Swiss Prot: Q00613

Purity: Affinity purification

For research use only.

Background:

The product of this gene is a heat-shock transcription factor. Transcription of heat-shock genes is rapidly induced after temperature stress. Hsp90, by itself and/or associated with multichaperone complexes, is a major repressor of this gene.

To place an order, please [Click HERE](#).