

Phospho-IKBKB-Y199

Reactivity: Human Mouse Rat

Tested applications: IF

Recommended Dilution: IF 1:100 - 1:200

Calculated MW: 87kDa

Observed MW: Refer to Figures

Immunogen:

A phospho specific peptide corresponding to residues surrounding Y199 of human IKBKB

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

ps

Synonym:

IKK2; IKKB; NFKB1KB; IKK-beta;

Catalog #: AP0369

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3551

Isotype: IgG

Swiss Prot: O14920

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene phosphorylates the inhibitor in the inhibitor/NF-kappa-B complex, causing dissociation of the inhibitor and activation of NF-kappa-B. The encoded protein itself is found in a complex of proteins. Several transcript variants, some protein-coding and some not, have been found for this gene.

To place an order, please [Click HERE](#).