

Phospho-JUNB-S79

Reactivity: Human Mouse Rat

Tested applications: IHC IF

Recommended Dilution: IHC 1:50 - 1:100 IF 1:100 - 1:200

Calculated MW: 43kDa

Observed MW: Refer to Figures

Immunogen:

A phospho specific peptide corresponding to residues surrounding S79 of human JUNB

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

k

Synonym:

AP-1;

To place an order, please [Click HERE](#).

Catalog #: AP0376

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3726

Isotype: IgG

Swiss Prot: P17275

Purity: Affinity purification

For research use only.