

Phospho-JUND-S255

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:1000 IHC 1:50 - 1:100

Calculated MW: 38kDa

Observed MW: Refer to Figures

Immunogen:

A phospho specific peptide corresponding to residues surrounding S255 of human JUND

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

k

Synonym:

AP-1

Catalog #: AP0377

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3727

Isotype: IgG

Swiss Prot: P17535

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this intronless gene is a member of the JUN family, and a functional component of the AP1 transcription factor complex. It has been proposed to protect cells from p53-dependent senescence and apoptosis. Alternate translation initiation site usage results in the production of different isoforms.

To place an order, please [Click HERE](#).