

## Phospho-Jun-T239

**Reactivity:** Human

**Tested applications:** WB

**Recommended Dilution:** WB 1:500 - 1:2000

**Calculated MW:** 43kDa

**Observed MW:** Refer to Figures

**Immunogen:**

A phospho specific peptide corresponding to residues surrounding T239 of human Jun

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

k

**Synonym:**

AP1; AP-1; c-Jun;

**Catalog #:** AP0378

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 3725

**Isotype:** IgG

**Swiss Prot:** P05412

**Purity:** Affinity purification

For research use only.

**Background:**

This gene is the putative transforming gene of avian sarcoma virus 17. It encodes a protein which is highly similar to the viral protein, and which interacts directly with specific target DNA sequences to regulate gene expression. This gene is intronless and is mapped to 1p32-p31, a chromosomal region involved in both translocations and deletions in human malignancies.

*To place an order, please [Click HERE](#).*