

Phospho-Jun-Y170

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 43kDa

Observed MW: Refer to Figures

Immunogen:

A phospho specific peptide corresponding to residues surrounding Y170 of human Jun

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

k

Synonym:

AP1; AP-1; c-Jun;

Catalog #: AP0381

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3725

Isotype: IgG

Swiss Prot: P05412

Purity: Affinity purification

For research use only.

Background:

This gene is the putative transforming gene of avian sarcoma virus 17. It encodes a protein which is highly similar to the viral protein, and which interacts directly with specific target DNA sequences to regulate gene expression. This gene is intronless and is mapped to 1p32-p31, a chromosomal region involved in both translocations and deletions in human malignancies.

To place an order, please [Click HERE](#).