

Phospho-MAP2K6-S207

Reactivity: Human Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:1000 IHC 1:50 - 1:100 IF 1:50 - 1:200

Calculated MW: 41kDa

Observed MW: Refer to Figures

Immunogen:

A phospho specific peptide corresponding to residues surrounding S207 of human MAP2K6

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

100µg/ml

Synonym:

MAP kinase kinase 6; MAP2K6; MAPK/ERK kinase 6; MAPKK 6; MEK6;

Background:

Catalyzes the concomitant phosphorylation of a threonine and a tyrosine residue in MAP kinase p38 exclusively.

To place an order, please [Click HERE](#).

Catalog #: AP0393

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5608

Isotype: IgG

Swiss Prot: P52564

Purity: Affinity purification

For research use only.