

Phospho-MARCKS-S162

Reactivity: Human Mouse Rat

Tested applications: WB IF

Recommended Dilution: WB 1:500 - 1:2000 IF 1:100 - 1:200

Calculated MW: 80kDa

Observed MW: Refer to Figures

Immunogen:

A phospho specific peptide corresponding to residues surrounding S162 of human MARCKS

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MACS; 80K-L; PKCSL; PRKCSL;

Catalog #: AP0403

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4082

Isotype: IgG

Swiss Prot: P29966

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a substrate for protein kinase C. It is localized to the plasma membrane and is an actin filament crosslinking protein. Phosphorylation by protein kinase C or binding to calcium-calmodulin inhibits its association with actin and with the plasma membrane, leading to its presence in the cytoplasm. The protein is thought to be involved in cell motility, phagocytosis, membrane trafficking and mitogenesis.

To place an order, please [Click HERE](#).