

Phospho-MET-Y1234/Y1235

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 155kDa

Observed MW: Refer to Figures

Immunogen:

A phospho specific peptide corresponding to residues surrounding Y1234/Y1235 of human MET

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HGFR; AUTS9; RCCP2; c-Met;

Catalog #: AP0408

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4233

Isotype: IgG

Swiss Prot: P08581

Purity: Affinity purification

For research use only.

Background:

The proto-oncogene MET product is the hepatocyte growth factor receptor and encodes tyrosine-kinase activity. The primary single chain precursor protein is post-translationally cleaved to produce the alpha and beta subunits, which are disulfide linked to form the mature receptor. Various mutations in the MET gene are associated with papillary renal carcinoma. Two transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).