

Phospho-PDPK1-S241

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100 IF 1:100 - 1:200

Calculated MW: 63kDa

Observed MW: Refer to Figures

Immunogen:

A phospho specific peptide corresponding to residues surrounding S241 of human PDPK1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PDK1; PDPK2; PRO0461;

To place an order, please [Click HERE](#).

Catalog #: AP0426

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5170

Isotype: IgG

Swiss Prot: O15530

Purity: Affinity purification

For research use only.