

Phospho-PIK3R1/PIK3R2/PIK3R3-Y467

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:1000

Calculated MW: 85kDa

Observed MW: Refer to Figures

Immunogen:

A phospho specific peptide corresponding to residues surrounding Y467 of human PIK3R1/2/3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PI3-kinase p85 subunit alpha ; PtdIns-3-kinase p85-alpha;

Catalog #: AP0427

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5295/ 8503/ 5296

Isotype: IgG

Swiss Prot: P27986/Q92569/O00459

Purity: Affinity purification

For research use only.

Background:

Binds to activated (phosphorylated) protein-Tyr kinases, through its SH2 domain, and acts as an adapter, mediating the association of the p110 catalytic unit to the plasma membrane. Necessary for the insulin-stimulated increase in glucose uptake and glycogen synthesis in insulin-sensitive tissues. Binds to activated (phosphorylated) protein-tyrosine kinases through its SH2 domain and regulates their kinase activity. During insulin stimulation, it also binds to IRS-1./Binds to activated (phosphorylated) protein-tyrosine kinases, through its SH2 domain, and acts as an adapter, mediating the association of the p110 catalytic unit to the plasma membrane.

To place an order, please [Click HERE](#).