

Phospho-PPP1R1B-T75

Reactivity: Human Mouse Rat

Tested applications: IF

Recommended Dilution: IF 1:100 - 1:200

Calculated MW: 32kDa

Observed MW: Refer to Figures

Immunogen:

A phospho specific peptide corresponding to residues surrounding T75 of human PPP1R1B

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DARPP32; DARPP-32;

Catalog #: AP0430

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 84152

Isotype: IgG

Swiss Prot: Q9UD71

Purity: Affinity purification

For research use only.

Background:

This gene encodes a bifunctional signal transduction molecule. Dopaminergic and glutamatergic receptor stimulation regulates its phosphorylation and function as a kinase or phosphatase inhibitor. As a target for dopamine, this gene may serve as a therapeutic target for neurologic and psychiatric disorders. Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).