

Phospho-PTEN-S380

Reactivity: Human Mouse Rat

Tested applications: IHC IF

Recommended Dilution: IHC 1:50 - 1:200 IF 1:100 - 1:200

Calculated MW: 47kDa

Observed MW: Refer to Figures

Immunogen:

A phospho specific peptide corresponding to residues surrounding S380 of human PTEN

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

ab

Synonym:

BZS; DEC; CWS1; GLM2; MHAM; TEP1; MMAC1; PTEN1; 10q23del;

Catalog #: AP0436

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5728

Isotype: IgG

Swiss Prot: P60484

Purity: Affinity purification

For research use only.

Background:

This gene was identified as a tumor suppressor that is mutated in a large number of cancers at high frequency. The protein encoded this gene is a phosphatidylinositol-3,4,5-trisphosphate 3-phosphatase. It contains a tensin like domain as well as a catalytic domain similar to that of the dual specificity protein tyrosine phosphatases. Unlike most of the protein tyrosine phosphatases, this protein preferentially dephosphorylates phosphoinositide substrates. It negatively regulates intracellular levels of phosphatidylinositol-3,4,5-trisphosphate in cells and functions as a tumor suppressor by negatively regulating AKT/PKB signaling pathway.

To place an order, please [Click HERE](#).