

Phospho-PXN-Y88

Reactivity: Human

Tested applications: WB IF

Recommended Dilution: WB 1:500 - 1:2000 IF 1:100 - 1:200

Calculated MW: 68kDa

Observed MW: Refer to Figures

Immunogen:

A phospho specific peptide corresponding to residues surrounding Y88 of human PXN

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

q

Catalog #: AP0441

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5829

Isotype: IgG

Swiss Prot: P49023

Purity: Affinity purification

For research use only.

Background:

This gene encodes a cytoskeletal protein involved in actin-membrane attachment at sites of cell adhesion to the extracellular matrix (focal adhesion). Alternatively spliced transcript variants encoding different isoforms have been described for this gene. These isoforms exhibit different expression pattern, and have different biochemical, as well as physiological properties (PMID:9054445).

To place an order, please [Click HERE](#).