

Phospho-PDK1-S241

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:63kDa

Observed MW:Refer to figures

Immunogen:

A phospho specific peptide corresponding to residues surrounding S241 of human PDK1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PDK1; PDPK2; PDPK2P; PRO0461;

Background:

Phosphoinositide-dependent kinase 1 (PDPK1, PDK1) is a serine/threonine protein kinase integral to the function of the PI 3-K/Akt signaling pathway. PDK1 and mTORC2 both phosphorylate and activate PKB/Akt, ensuring a cellular response to stimuli such as growth factors and insulin signaling. Akt is the main effector of PDK1.

To place an order, please [Click HERE](#).

Catalog #:AP0477

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:5170

Isotype:IgG

Swiss Prot:O15530

Purity:Affinity purification

For research use only.