

Phospho-KSR1-S392

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:102kDa

Observed MW:Refer to figures

Immunogen:

A phospho specific peptide corresponding to residues surrounding S392 of human KSR1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

KSR; RSU2;

To place an order, please [Click HERE](#).

Catalog #:AP0487

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:8844

Isotype:IgG

Swiss Prot:Q81VT5

Purity:Affinity purification

For research use only.