

Phospho-KSR1-S404

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 102kDa

Observed MW: Refer to figures

Immunogen:

A phospho specific peptide corresponding to residues surrounding S404 of human KSR1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

KSR; RSU2;

To place an order, please [Click HERE](#).

Catalog #: AP0488

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8844

Isotype: IgG

Swiss Prot: Q8IVT5

Purity: Affinity purification

For research use only.