

b NGF Human, HEK

Description: Beta-NGF Human Recombinant produced in HEK cells is a non-glycosylated non-disulfide linked homodimer, having a total molecular weight of 13kDa. The b-NGF is purified by proprietary chromatographic techniques.

Catalog #: CYP5-086

For research use only.

Synonyms: Beta Polypeptide, NGF, NGFB, HSN5, Beta-NGF, MGC161426, MGC161428.

Source: HEK.

Physical Appearance: Sterile Filtered White lyophilized (freeze-dried) powder.

Purity: Greater than 95% as observed by SDS-PAGE.

Formulation:

The b-NGF was lyophilized from 1mg/ml in 1xPBS.

Stability:

Lyophilized b-NGF although stable at room temperature for 3 weeks, should be stored desiccated below -18°C. Upon reconstitution b-NGF should be stored at 4°C between 2-7 days and for future use below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze-thaw cycles.

Usage:

NeoBiolabs products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Solubility:

It is recommended to reconstitute the lyophilized b-NGF in sterile water not less than 100

Introduction:

NGF-beta has nerve growth stimulating activity and the complex is involved in the regulation of growth and the differentiation of sympathetic and certain sensory neurons. Mutations in this gene have been associated with hereditary sensory and autonomic neuropathy, type 5 (HSAN5), and dysregulation of this gene's expression is associated with allergic rhinitis.

Biological Activity:

The specific activity was determined by the dose-dependent stimulation of the proliferation of human TF-1 cells (human erythroleukemic indicator cell line) and is typically 1.5-8ng/ml.

To place an order, please [Click HERE](#).