

IFNA7 Human

Description: IFNA7 Human Recombinant produced in E.Coli is a single, non-glycosylated polypeptide chain containing 191 amino acids (24-189 a.a) and having a molecular mass of 22.3kDa. IFNA7 is fused to a 25 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques.

Catalog #: CYP5-203

For research use only.

Synonyms: Interferon alpha-7, IFN-alpha-7, Interferon alpha-J, LeIF J, Interferon alpha-J1, IFN-alpha-J1, IFNA7, IFNA-J, IFN-alphaJ.

Source: E.coli.

Physical Appearance: Sterile Filtered colorless solution.

Amino Acid Sequence: MGSSHHHHHH SSGLVPRGSH MGSHMCDLPQ THSLRNRRAL
ILLAQMGRIS PFSCCLKDRHE FRFPEEEFDG HQFQKTQAIS VLHEMIQQTFLNFSTEDSSA
AWEQSLLEKF STELYQQLND LEACVIEVG VEETPLMNED FILAVRKYFQ RITLYLMEKK
YSPCAWEVVR AEIMRSFSFS TNLKKGLRRK D.

Purity: Greater than 90% as determined by SDS-PAGE.

Formulation:

IFNA7 protein solution (0.5mg/ml) containing 20mM Tris-HCl buffer (pH 8.0), 0.4M urea and 10% glycerol.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

Interferon alpha 7 (IFNA7) is a member of the alpha/beta interferon family. IFNA7 is generated by macrophages. IFN-alpha has antiviral functions. Interferon promotes the production of 2 enzymes: a protein kinase and an oligoadenylate synthetase.

To place an order, please [Click HERE](#).