

sIL 2R Human

Description: Interleukin-2R Human Recombinant produced in E.Coli is a single, non-glycosylated polypeptide chain containing 192 amino acids (22-213) corresponding to the mature IL2-R protein with an amino terminal hexahistidine tag. The IL-2R is purified by proprietary chromatographic techniques.

Catalog #: CYPs-292

For research use only.

Synonyms: CD25, IL2R, TCGFR, IL-2RA, sIL-2RA, TAC antigen, sIL-2R, IDDM10, p55.

Source: Escherichia Coli.

Physical Appearance: Sterile colorless liquid formulation.

Purity: Greater than 95.0% as determined by (a) Analysis by RP-HPLC. (b) Analysis by SDS-PAGE.

Formulation:

IL-2R His-Tag protein is supplied in 1x PBS and 50% glycerol pH 7.4.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Please avoid freeze thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

IL2-Ra is one of the three constituent subunits of the IL2 receptor. IL-2Ra is released into the serum after increased cellular expression such as increased activation of B and T cells. Clinical manifestations of IL2-Ra elevation include autoimmune conditions and some leukemias and lymphomas.

To place an order, please [Click HERE](#).