

Leptin qA Ovine

Description:Leptin Antagonist Quadruple Mutant Ovine Recombinant is a single non-glycosylated polypeptide chain containing 146 amino and additional Ala at N-terminus acids and having a molecular mass of ~ 16 kDa, Leptin was mutated, resulting in L39A/D40A/F41A/I42A mutant. Leptin Antagonist Quadruple Mutant Ovine Recombinant was purified by proprietary chromatographic techniques.

Catalog #:CYPS-364

For research use only.

Source:Escherichia coli.

Physical Appearance:White lyophilized (freeze-dried) powder.

Amino Acid Sequence:The sequence of the first five N-terminal amino acids was determined and was found to be Ala-Val-Pro-Ile-Arg.

Purity:Greater than 98.0% as determined by:(a) Gel filtration analysis.(b) Analysis by SDS-PAGE.

Formulation:

The protein was lyophilized from a concentrated (1mg/ml) solution with 0.0045mM NaHCO₃.

Stability:

Lyophilized Leptin-Antagonist Quadruple Mutant Ovine Recombinant although stable at room temperature for several weeks, should be stored desiccated below -18°C. Upon reconstitution at > 0.1 LEP-qA mutant mg/ml and up to 2 mM and filter sterilization Leptin mutant can be stored at 4°C or even room temperature for several weeks making it suitable for long term infusion studies using osmotic pumps. At lower concentration addition of a carrier protein (0.1% HSA or BSA) is suggested. Please prevent freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Solubility:

It is recommended to reconstitute the lyophilized Leptin-Antagonist Quadruple Mutant Ovine Recombinant in sterile 0.4% NaHCO₃ adjusted to pH 8-9, not less than 100

Biological Activity:

NeoBiolabs leptin qA is capable of inhibiting leptin-induced proliferation of BAF/3 cells stably transfected with the long form of ovine leptin receptor. It also inhibits various leptin effects in several in vitro bioassays.

To place an order, please [Click HERE](#).