

IL 7 Mouse

Description: Interleukin-7 Mouse Recombinant produced in E.Coli is a single, non-glycosylated polypeptide chain containing 130 amino acids and having a molecular mass of 15.0 kDa. The IL-7 is purified by proprietary chromatographic techniques.

Catalog #: CYPs-379

For research use only.

Synonyms: Lymphopoietin 1 (LP-1), pre-B cell factor, IL-7.

Source: Escherichia Coli.

Physical Appearance: Sterile Filtered White lyophilized (freeze-dried) powder.

Amino Acid Sequence: MECHIKDKEG KAYESVLMIS IDELDKMTGT DSNCPNNEPN
FFRKHVCDT KEAAFLNRRA RKLKQFLKMN ISEEFNVHLL TVSQGTQTLV NCTSKEEKNV
KEQKKNDACF LKRLLEIKT CWNKILKGS.

Purity: Greater than 95.0% as determined by: (a) Analysis by RP-HPLC. (b) Analysis by SDS-PAGE.

Formulation:

Filtered (0.2

Stability:

Lyophilized Interleukin-7 although stable at room temperature for 3 weeks, should be stored desiccated below -18°C. Upon reconstitution IL7 should be stored at 4°C between 2-7 days and for future use below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Solubility:

It is recommended to reconstitute the lyophilized Interleukin 7 in sterile 18M-cm H₂O not less than 100µg/ml, which can then be further diluted to other aqueous solutions.

Introduction:

Interleukin-7 is a potent lymphoid cell growth factor produced primarily by stromal cells. IL-7 has been shown to support the proliferation and differentiation of pre B- and early T cells as well as displaying a biological effect on cells of NK and myeloid lineages.

Biological Activity:

The ED₅₀ as determined by the dose-dependant stimulation of murine 2E8 cells is < 0.2ng/ml, corresponding to a Specific Activity of 5,000,000 units/mg.

To place an order, please [Click HERE](#).