

IL 1 beta Human, His

Description: Interleukin-1 beta His-Tag Human Recombinant produced in E.Coli is a single, non-glycosylated, Polypeptide chain containing 153 amino acids fragment (117-269) and having a total molecular mass of 21.88 kDa with an amino-terminal hexahistidine tag. The IL-1b His is purified by proprietary chromatographic techniques.

Catalog #: CYPs-487

For research use only.

Synonyms: Catabolin, Lymphocyte-activating factor (LAF), Endogenous Pyrogen (EP), Leukocyte Endogenous Mediator (LEM), Mononuclear Cell Factor (MCF), IL1F2, IL-1 beta.

Source: Escherichia Coli.

Physical Appearance: Sterile Filtered clear solution.

Purity: Greater than 95.0% as determined by (a) Analysis by RP-HPLC. (b) Analysis by SDS-PAGE.

Formulation:

IL-1b His is supplied in 1x PBS and 50% glycerol.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Please avoid freeze thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

Interleukin-1b is a potent pro-inflammatory cytokine produced by a wide variety of cell types including monocytes and macrophages. It displays a broad range of biological activity including activation of B and T cells in response to inflammation and the activation of endothelial cells.

To place an order, please [Click HERE](#).