

Leptin Receptor Chicken

Description:Leptin Binding Domain Chicken Recombinant also called Leptin Receptor produced in E.Coli is a single, non-glycosylated, polypeptide chain containing 208 amino acids and having a molecular mass of 24.5 kDa. Chicken Leptin Receptor consists of the cytokine binding domain of leptin receptor amino acids 420-626 of chicken leptin receptor. The Leptin Binding Domain is purified by proprietary chromatographic techniques.

Catalog #:CYP5-516

For research use only.

Synonyms:OB Protein, Obesity Protein, OBS, Obesity factor, Leptin Receptor.

Source:Escherichia Coli.

Physical Appearance:Sterile Filtered colorless solution at a concentration of 0.4 mg/ml.

Amino Acid Sequence:The sequence of the first six N-terminal amino acids was determined and was found to be Ala-Ile-Asp-Val-Asn-Ile. Biological Activity is evidenced by high affinity binding of mammalian leptins at 1:1 molar ratio.

Purity:Greater than 99.0% as determined by:(a) Analysis by SEC-HPLC.(b) Analysis by SDS-PAGE.

Formulation:

The protein was filter sterilized and stored at 4°C (0.2 to 0.5 mg/ml) solution of Tris-HCl buffer, pH 9.0 with 150mM NaCl.

Stability:

Sterile solutions at 0.5mg/ml or less are stable at 4°C for several months.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

To place an order, please [Click HERE](#).