

SHH Human, His

Description: Sonic HedgeHog Recombinant Human produced in E.Coli is a single, non-glycosylated polypeptide chain containing 183 amino acids and having a molecular mass of 20.7kDa. The Sonic HedgeHog is fused to 8 amino acid His-Tag at C-terminus and purified by proprietary chromatographic techniques.

Catalog #: CYPs-724

For research use only.

Synonyms: SHH, HHG-1, HHG1, Sonic hedgehog protein, TPT, HLP3, HPE3, SMMCI, TPTPS, MCOPCB5.

Source: Escherichia Coli.

Physical Appearance: Sterile filtered colorless solution.

Amino Acid Sequence: MCGPGRGFGK RRHPKKLTPL AYKQFIPNVA EKT LGASGRY
EGKISRNSER FKELTPNYP DIIFKDEENT GADRLMTQRC KDKLNALAI SVMNQWPGVKL
RVTEGWDEDG HHSEESLHYE GRAVDITSD RDRSKYGMLA RLAVEAGFDW VYYESKAHIH
CSVKAENSVA AKSGGLEHHH HHH.

Purity: Greater than 95.0% as determined by SDS-PAGE.

Formulation:

SHH Human solution containing 20mM Tris HCL pH-8, 0.1M NaCl, & 10% glycerol.

Stability:

SHH Human although stable at 4°C for 1 week, should be stored desiccated below -18°C. Please prevent freeze thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

Recombinant Human Sonic Hedgehog is part of a small group of secreted proteins that are vital for development in both vertebrates and invertebrates. 3 mammalian hedgehog genes (sonic, desert, Indian) share about 60% homology. The Human Sonic Hedgehog is 99% homologous to the mouse gene. Sonic HedgeHog is a protein that is vital in guiding the early embryo. It has been associated as the major inductive signal in patterning of the ventral neural tube, the anterior-posterior limb axis, and the ventral somites. Sonic HedgeHog binds to the patched receptor, which functions in association with smoothened, to activate the transcription of target genes. In the absence of sonic HedgeHog, patched receptor represses the constitutive signaling activity of smoothened. Sonic HedgeHog also regulates another factor, the gli oncogene. Sonic HedgeHog intercellular signal is essential for a various patterning events during development: signal produced by the notochord that induces ventral cell fate in the

To place an order, please [Click HERE](#).