

Dengue NS1 c

Description: The E.Coli derived recombinant protein contains the NS1 Dengue Virus c-end Type-2 immunodominant regions.

Catalog #: DEPS-009

Purity: Dengue Protein is >90% pure as determined by 10% PAGE (coomassie staining).

For research use only.

Purification Method:

Dengue protein is purified by proprietary chromatographic technique.

Formulation:

20mM Tris-HCl pH 8.0, 8M urea & 60mM NaCl.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. They may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Applications:

Each laboratory should determine an optimum working titer for use in its particular application.

Introduction:

Caused by one of four closely related virus serotypes of the genus Flavivirus, family Flaviviridae, each serotype is sufficiently different that there is no cross-protection and epidemics caused by multiple serotypes (hyperendemicity) can occur. In cell culture experiments and mice Morpholino antisense oligos have shown specific activity against Dengue virus.

Storage:

Dengue although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.

To place an order, please [Click HERE](#).