

Dengue NS1

Description: The E.Coli derived recombinant protein contains the NS1 Dengue Virus full length Type-2 immunodominant regions, having an Mw of 45kDa. The dengue protein is fused to 6xHis tag and purified by proprietary chromatographic technique.

Catalog #:DEPS-011

For research use only.

Purity: Protein is >95% pure as determined by SDS-PAGE.

Formulation:

Phosphate buffered saline, pH 7.4, 100mM arginine and 2mM DTT.

Stability:

Dengue NS1 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. They may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Applications:

Each laboratory should determine an optimum working titer for use in its particular application.

Introduction:

Caused by one of four closely related virus serotypes of the genus Flavivirus, family Flaviviridae, each serotype is sufficiently different that there is no cross-protection and epidemics caused by multiple serotypes (hyperendemicity) can occur. In cell culture experiments and mice Morpholino antisense oligos have shown specific activity against Dengue virus.

To place an order, please [Click HERE](#).