

Dengue Envelope-3 22kDa

Description: The E.coli derived recombinant 28 kDa protein is genetically engineered peptide which is derived from Dengue Type-3 Envelope. This region also contains a common antigen for Dengue type 1, 2 and 4. The sequence is from Dengue virus type 3 strain C0360/94, genebank access code is AY923865.

Catalog #:DEPS-015

For research use only.

Purity: Protein is >95% pure as determined by 12% PAGE (coomassie staining).

Formulation:

Phosphate buffered saline, pH-7.4.

Stability:

Dengue Envelope ST3 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. They may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Applications:

Each laboratory should determine an optimum working titer for use in its particular application.

Introduction:

Caused by one of four closely related virus serotypes of the genus Flavivirus, family Flaviviridae, each serotype is sufficiently different that there is no cross-protection and epidemics caused by multiple serotypes (hyperendemicity) can occur. In cell culture experiments and mice Morpholino antisense oligos have shown specific activity against Dengue virus.

To place an order, please [Click HERE](#).