

## Dengue NS1 ST1, Insect

**Description:** Recombinant Dengue Virus NS1 Subtype 1 produced in Insect Cells is a polypeptide chain containing amino acids 777-1131 and having a molecular weight of approximately 50kDa. Dengue NS1 ST1 is purified by proprietary chromatographic technique.

Catalog #: DEPS-036

For research use only.

**Source:** Insect cells.

**Physical Appearance:** Sterile filtered colorless solution.

**Purity:** Protein is >95% pure as determined by 12.5% SDS-PAGE.

**Formulation:**

Dengue NS1 ST1 protein solution (1mg/ml) in 1xPBS, pH7.4, 0.1% Thimerosal, 5mM EDTA, 1

**Stability:**

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

**Usage:**

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. They may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

**Applications:**

Ag test control.

**Introduction:**

Caused by one of four closely related virus serotypes of the genus *Flavivirus*, family *Flaviviridae*, each serotype is sufficiently different that there is no cross-protection and epidemics caused by multiple serotypes (hyperendemicity) can occur. In cell culture experiments and mice Morpholino antisense oligos have shown specific activity against Dengue virus.

**To place an order, please [Click HERE](#).**