

LCN2 Human, Pichia

Description: Neutrophil Gelatinase Associated Lipocalin Human Recombinant is expressed in Pichia pastoris having a molecular weight of approximately 21kDa. The LCN2 is purified by proprietary chromatographic techniques.

Catalog #: ENPS-071

For research use only.

Synonyms: Neutrophil gelatinase-associated lipocalin, NGAL, p25, 25 kDa
alpha-2-microglobulin-related subunit of MMP-9, Lipocalin-2, Oncogene 24p3, LCN2.

Source: Pichia pastoris.

Physical Appearance: Sterile Filtered White lyophilized (freeze-dried) powder.

Purity: Greater than 96% as determined by SDS-PAGE.

Formulation:

LCN2 was lyophilized from 0.02M NH₄HCO₃.

Stability:

Lyophilized LCN2 although stable at room temperature for 3 weeks, should be stored desiccated below -18°C. Upon reconstitution LCN2 should be stored at 4°C between 2-7 days and for future use below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Solubility:

It is recommended to reconstitute the lyophilized LCN2 in sterile 18M-cm H₂O not less than 100µg/ml, which can then be further diluted to other aqueous solutions.

Introduction:

Recombinant Human Lipocalin-2 also called Neutrophil Gelatinase Associated Lipocalin (NGAL) belongs to a family of lipocans which include 25 proteins (including α1-microglobulin and b-lactoglobulin), which are characterized by their ability to bind small lipophilic substances in their hydrophobic core. They thereby serve as transporters of substances like retinal, biliverdins & prostaglandins. There are indications that NGAL is involved in modulation of the inflammatory response and is found in the plasma of patients after stroke.

To place an order, please [Click HERE](#).