

MSRB3 Human

Description:MSRB3 Human Recombinant fused with an 8 amino acid His tag at C-terminus produced in E.Coli is a single, non-glycosylated, polypeptide chain containing 174 amino acids (21-185 a.a.) and having a molecular mass of 19kDa. The MSRB3 is purified by proprietary chromatographic techniques.

Catalog #:ENPS-100

For research use only.

Synonyms:Methionine-R-sulfoxide reductase B3, MSRB3, DFNB74, FLJ36866, DKFZp686C1178.

Source:Escherichia Coli.

Physical Appearance:Sterile Filtered colorless solution.

Amino Acid Sequence:MCGLPSGSCR DKKNCKVVFS QQELRKRLTP LQYHVTQKEG
TESAFEGEYT HHKDPGIYKC VVCGTPLFKS ETKFDSGSGW PSFHDVINSE AITFTDDFSY
GMHRVETSCS QCGAHLGHIF DDGPRPTGKR YCINSAALSF TPADSSGTAE GGSGVASPAQ
ADKAELLEHH HHHH.

Purity:Greater than 90.0% as determined by SDS-PAGE.

Formulation:

The MSRB3 solution (1 mg/ml) contains 20mM Tris-HCl buffer (pH 8.0), 1mM DTT, 10% glycerol and 0.1mM PMSF.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

Methionine sulfoxide reductase B3 (MSRB3) is a member of the methionine sulfoxide reductases (MSR) family proteins. MSRB3 catalyzes the reduction of methionine sulfoxide to methionine. The MSRB3 enzyme acts as a monomer and requires zinc as a cofactor. MSRs are thought to defend against reactive oxygen species-induced oxidative damage in various organs, including the most environmentally exposed organ, the human skin. MSRB3 has a vital role in cold tolerance by eliminating MetO and ROS which accumulate at the ER during cold acclimation.

To place an order, please [Click HERE](#).