

TREX2 Human

Description: TREX2 produced in E.Coli is a single, non-glycosylated polypeptide chain containing 256 amino acids (1-236a.a.) and having a molecular mass of 28.0 kDa. TREX2 is fused to a 20 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques.

Catalog #: ENPS-102

For research use only.

Synonyms: Three Prime Repair Exonuclease 2, 3'-5' exonuclease TREX2 long form.

Source: Escherichia Coli.

Physical Appearance: Sterile Filtered clear solution.

Amino Acid Sequence: MGSSHHHHH SSGLVPRGSH MSEAPRAETF VFLDLEATGL
PSVEPEIAEL SLFAVHRSSL ENPEHDESGA LVLPRVLDKL TLCMCPEPFP TAKASEITGL
SSEGLARCRK AGFDGAVVRT LQAFLSRQAG PICLVAHNGF DYDFPLLCAE LRRLGARLPR
DTVCLDTLPA LRGLDRAHSH GTRARGRQGY SLGSLFHRYF RAEPSAAHSA EGDVHTLLLI
FLHRAAELLA WA

Purity: Greater than 95% as determined by SDS-PAGE.

Formulation:

The TREX2 protein solution (1mg/1ml) is formulated in 20mM Tris-HCl buffer (pH8.0), 200mM NaCl, 5mM DTT and 30% glycerol.

Usage:

NeoBiolabs products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

TREX2 holds a 3-prime-to-5-prime exonuclease activity and eliminates mismatched, modified, fragmented, and normal nucleotides to produce the appropriate 3-prime termini for following steps in the DNA metabolic pathways. TREX2 has a role in DNA replication, repair, and recombination.

Storage:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

To place an order, please [Click HERE](#).