

NUDT10 Human

Description:NUDT10 produced in E.Coli is a single, non-glycosylated polypeptide chain containing 172 amino acids (1-164 a.a.) and having a molecular mass of 19.5kDa.NUDT10 is fused to an 8 amino acid His-tag at C-terminus & purified by proprietary chromatographic techniques.

Catalog #:ENPS-133

For research use only.

Synonyms:Diphosphoinositol polyphosphate phosphohydrolase 3-alpha, DIPPP-3-alpha, DIPPP3-alpha, hDIPPP3alpha, Diadenosine 5',5'''-P1,P6-hexaphosphate hydrolase 3-alpha, Nucleoside diphosphate-linked moiety X motif 10, Nudix motif 10, hAps2, NUDT10, APS2, DIPPP3A.

Source:Escherichia Coli.

Physical Appearance:Sterile Filtered colorless solution.

Amino Acid Sequence:MKCKPNQTRT YDPEGFKKRA ACLCFRSERE DEVLVSSSR
YPDRWIVPGG GMEPEEEPPG AAVREYEEA GVKGKLGRL GVFEQNQDPK HRTYVYVLT
TELLEDWEDS VSIGRKREWF KVEDAIKVLQ CHKPVHAEYL EKLKLGGSPT NGNSMAPSSP
DSDPLEHHHH HH.

Purity:Greater than 95.0% as determined by SDS-PAGE.

Formulation:

NUDT10 protein solution (0.5mg/ml) containing 20mM Tris-HCl buffer (pH8.0), 20% glycerol, 0.1M NaCl and 1mM DTT.

Stability:

NUDT10 Human Recombinant although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

NUDT10 belongs to the Nudix hydrolase family of pyrophosphatases. Nudix hydrolases contain a characteristic Nudix domain and are responsible for catalyzing the hydrolysis of nucleoside diphosphate derivatives. NUDT10 functions as a manganese-dependent polyphosphate phosphohydrolase with an optimum pH of 8.5. NUDT10 protein specifically metabolizes diadenosine-polyphosphates and, to a lesser extent, diphosphoinositol polyphosphates.

To place an order, please [Click HERE](#).