

SEN8 Human

Description: SEN8 Human Recombinant produced in E.coli is a single, non-glycosylated polypeptide chain containing 232 amino acids (1-212) and having a molecular mass of 26.3 kDa. The SEN8 is fused to a 20 amino acid His-Tag at N-terminus and purified by proprietary chromatographic techniques.

Catalog #: ENPS-153

For research use only.

Synonyms: SUMO/sentrin specific peptidase family member 8, DEN1, NEDP1, Protease cysteine 2 (NEDD8 specific), PRSC2, NEDD8-specific protease 1, HsT17512, Deneddylase-1, NEDD8 specific-protease cysteine 2, Sentrin/SUMO-specific protease SEN8.

Source: Escherichia Coli.

Physical Appearance: SEN8 is supplied as a sterile filtered clear solution.

Amino Acid Sequence: MGSSHHHHHH SSGLVPRGSH MDPVLSYMD SLLRQSDVSL
LDPPSWLNDH IIGFAFEYFA NSQFHDCSDH VSFISPEVTQ FIKCTSNPAE IAMFLEPLDL
PNKRVVFLAI NDNSNQAAGG THWSLLVYLQ DKNSFFHYDS HSRSNSVHAK QVAEKLEAFL
GRKGDKLAFV EEKAPAQONS YDCGMYVICN TEALCQNFFR QQTESLLQLL TPAYITKKRG
EWKDLITTLA KK

Purity: Greater than 95.0% as determined by SDS-PAGE.

Formulation:

SEN8 protein 1mg/ml is supplied in 20mM Tris-HCL, pH-8, 0.1M NaCl, 1mM DTT and 10% Glycerol.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

SEN8 is a cysteine protease which belongs to the sentrin-specific protease family. SEN8 takes part in processing and deconjugation of the ubiquitin-like protein labeled, neural precursor cell expressed developmentally downregulated 8 (NEDD8).

To place an order, please [Click HERE](#).