

LCAT Human

Description: LCAT Human Recombinant produced in E.Coli is a single, non-glycosylated, polypeptide chain containing 441 amino acids (25-440) which includes a 20 amino acid His Tag fused at N-terminus and having a total molecular mass of 49.8 kDa. LCAT Human Recombinant is purified by proprietary chromatographic techniques.

Catalog #: ENPS-387

For research use only.

Synonyms: Phosphatidylcholine-sterol acyltransferase, Lecithin-cholesterol acyltransferase, Phospholipid-cholesterol acyltransferase, LCAT.

Source: Escherichia Coli.

Physical Appearance: Sterile filtered colorless solution.

Amino Acid Sequence: MGSSHHHHHH SSGLVPRGSH MGSHMFLLN VLFPPHTTPK
AELSNHTRPV ILVPGCLGNQ LEAKLDKPDV VNWMCYRKTE DFFTIWLDLN MFLPLGVDCW
IDNTRVYVNR SSGLVSNAPG VQIRVPGFGK TYSVEYLDSS KLAGYLHTLV QNLVNNGYVR
DETVRAAPYD WRLEPGQQEEYRKLGLVE EMHAAYGKPV FLIGHSLGCL HLLYFLLRQP
QAWKDRFIDG FIS

Purity: Greater than 85.0% as determined by SDS-PAGE.

Formulation:

The LCAT 0.5mg/ml solution contains 20mM Tris-HCl buffer (pH 8.0), 0.5mM DTT, 0.4M Urea and 10% glycerol.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Usage:

NeoBiolab's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Introduction:

LCAT is an extracellular cholesterol esterifying enzyme, lecithin-cholesterol acyltransferase. The esterification of cholesterol is required for cholesterol transport. LCAT is an essential enzyme in the extracellular metabolism of plasma lipoproteins.

To place an order, please [Click HERE](#).